



TR Scan

Anti-vibration system

TR Scan

- ✓ Basic solution
- ✓ Passive air solution
- ✓ Dedicated active solution
- ✓ Custom solution




DESCRIPTION

The main enemies when performing a measurement are external vibrations.
Trimos offers 3 solutions to minimise the influence of vibrations.

Antivibration pad

TA-TO-604



Basic solution :
Rubber plate with 3 layers of absorption.
Surface delivered 200x200 mm to be cut.
Very easy to install under TR SCAN feet.

Air passive pad

TA-TO-603



Passive solution FAEBI50
Replaces the 3 standard feet.
Bears TR SCAN weight with parts up to 25 kg.
Elevation 60 mm.
Base width 110 mm.

Dedicated Active antivibration table


TA-TO-601



Active solution dedicated to TR SCAN :
Very rigid table equipped with wheels, as well as with a antivibration plate dedicated to TR SCAN, with working surface.
Length : 1300 mm
Width : 900 mm
Height : 750 mm

Solution

On request



Custom solution on request.

